



AESCULAP®

All it takes to operate.

B | BRAUN
SHARING EXPERTISE

AESCULAP AG & CO. KG

Am Aesculap-Platz
78532 Tuttlingen

Phone +49 7461 95-0

Fax +49 7461 95-2600


www.aesculap.de

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.

B. Braun Sutures

Surgical Loop

For retraction of organs



Special Suture Set made of silicone, cotton or polyester

Surgical Loop

For retraction of organs

Surgical Loop are retraction tapes made of radiopaque silicone, cotton or polyester. The surgical tapes are used for the intraoperative isolation, marking, mechanical support and looping of organs, blood vessels, tendons and nerves, or for tying off the umbilical cord of the new born infant.

Surgical retraction tape made of silicone, cotton or polyester

- To identify organs/areas while operating and to retract or hold nerves, tendons, arteries, veins, ureters and vessels in order to ease the access to the operation site
- To clamp
- To tie off the umbilical cord of the newborn infant

Surgical Loop ensures and provides:

- Flexibility and resistance
- Smooth surface with no tissue adherence
- Attractive colour presentation to ease differentiation
- Different product dimensions to suit all purposes
- Silicone radiopaque material to ensure patient safety



silicone, cotton or polyester

non-absorbable

uncoated

tube or braided tape

Product range

Radiopaque silicone tube

Colour / Code

Diameter	Length	Units / box	White	Blue	Yellow	Red
1.5 mm	2x45 cm	24	B1095048	B1095021	B1095030	B1095013
1.5 mm	75 cm	24	B1095242	B1095226	B1095234	B1095218
2.5 mm	2x45 cm	24	B1095145	B1095129	B1095137	B1095110
2.5 mm	75 cm	24	B1095544	B1095528	B1095536	B1095510
1.2 mm	2x45 cm	24	B1095552	The thread is attached to a DRN 25 needle		DRN 25

Cotton braided tape

Colour / Code

Width	Length	Units / box	Purple	Blue	White	Green
4 mm	2x38 cm	24			B1094920	
4 mm	75 cm	24	B1095641	B1095633	B1095625	B1095650
8 mm	75 cm	24		B1095668	B1095676	
4 mm	75 cm	24	B1095617	The tape is attached to a GS 65 needle		GS 65

Polyester braided tape

Colour / Code

Width	Length	Units / box	Green
3.2 mm	2x75 cm	24	B1094904
4 mm	75 cm	24	B1094912

Product features

Composition	Silicone, cotton or polyester
Structure	Tube or braided tape
Coating	Uncoated
Color	Red, blue, yellow, white, green or violet
Sizes	Tube: 1.2 mm, 1.5 mm or 2.5 mm Tape: 3.2 mm, 4 mm or 8 mm
Mass absorption	Non-absorbable
Indications	Organ retraction
Contraindications	Not known
Sterilisation	Ethylene oxide

